BRA~~X~~ON<sup>®</sup>  
*FAST*

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
---

## **BRAXON® fast ACELLULAR DERMAL MATRIX FOR BREAST RECONSTRUCTION**

### **DESCRIPTION OF THE DEVICE**

BRAXON® fast is a non-pyrogenic pre-shaped porcine dermis-derived collagen membrane (acellular dermal matrix - ADM). This medical device is resistant and biocompatible and is embedded in the patient's tissue through cellular and microvascular infiltration.

### **INDICATIONS**

The BRAXON® fast membrane is specifically designed for prepectoral implant in breast surgery and, in particular, post-mastectomy breast reconstruction, revision surgery, and cosmetic surgery.

BRAXON® fast has a patented dome-shaped three-dimensional shape, which allows for housing a round or anatomical breast implant of up to 600 g. This membrane allows for subcutaneous placement of the breast implant without having to move the pectoralis major muscle, thereby maintaining its functionality.

### **CONTRAINDICATIONS**

BRAXON® fast is not designed, sold, or intended to be used in any way other than as described in the indications. It is not to be used in the following circumstances:

- patients with known hypersensitivity to porcine collagen
- acute, chronic infections
- autoimmune soft tissue or connective tissue diseases
- implantation in ischemic or necrotic areas or in patients with altered microcirculation.

### **WARNINGS AND STORAGE INSTRUCTIONS**

- Keep away from heat sources or direct sunlight.
- Store the product in its original packaging in a dry and clean place at room temperature.
- Do not use this product after its expiration date.
- This is a disposable product. Do not sterilise.
- Do not use in the event of open or damaged packaging.
- Do not use if the material or the sterile internal pouch have come into contact with non-sterile surfaces or objects.
- After use, dispose of all the unused parts of the product and its packaging materials in compliance with recognised medical practices and the national and regional environmental laws in force concerning the disposal of biological material and packaging.
- This device must be used solely by specialist medical personnel duly trained in tissue regeneration techniques.

### **STERILISATION**

This product is sterilised with ethylene oxide (EtO) and cannot be re-sterilised. Sterility is guaranteed

only if the product is in its intact and undamaged original packaging. Use during pregnancy and breastfeeding

There are currently no studies on the use of the BRAXON® fast membrane during pregnancy or breastfeeding or its impact on human fertility. Therefore, the benefits for the mother and possible risks for the baby must be assessed case by case before using the BRAXON® fast matrix.

### **USE DURING PREGNANCY AND BREASTFEEDING**

There are currently no studies on the use of the BRAXON® fast membrane during pregnancy or breastfeeding or its impact on human fertility. Therefore, the benefits for the mother and possible risks for the baby must be assessed case by case before using the BRAXON® fast matrix.

### **SIDE EFFECTS**

As with all medical devices, there may be possible complications along with the benefits. Potential side effects typically associated with acellular matrices are the same as those reported in breast surgery, including redness, inflammation, serous secretions, infection, formation of haematomas and seromas, and recurrent tissue damage.

### **INFORMATION**

The BRAXON® fast membrane is of natural origin; therefore, when dry,

there may be some rippling in the collagen structure and a variation in the membrane's colour and thickness. These phenomena do not affect the quality and function of BRAXON® fast.

The BRAXON® fast matrix does not contain preservatives and is supplied dehydrated in a sterile packaging.

The BRAXON® fast membrane is available in three sizes (small, medium, and large) suitable for containing implants with different width, projection, and volume (see size chart).

The BRAXON® fast membrane can be cut out and adapted to the size of the breast implant both in its dehydrated and hydrated state using a pair of sterile scissors.

In the event of contamination or severe infection, it is essential to take the measures provided for by institutional infection control protocols. (The presence of a significant bacterial load may compromise the functionality of the surgical membrane. The use of antimicrobial drugs could slow down the tissue regeneration process).

### **INSTRUCTIONS FOR IMPLANTATION**

We recommend performing a pinch test of at least 1 cm before the surgery to ensure sufficient thickness of the subcutaneous layer. We recommend assessing the vascularisation and vitality of the mastectomy flap even after removing the breast. The acellular

matrix should be placed in close contact with healthy and well-vascularised tissue to promote cell migration and tissue regeneration.

1. Open the box and remove the double sterile packaging.
2. Open the external wrapping and remove the internal pouch under sterile conditions. The internal pouch is sterile and can be placed directly in the sterile area of the surgery.
3. Carefully open the internal pouch and remove BRAXON® fast under sterile conditions.
4. Handle BRAXON® fast always using sterile gloves and tweezers.
5. BRAXON® fast is supplied sterile and must be rehydrated before implantation to facilitate suturing. Place BRAXON® fast in a sterile bowl filled with room temperature saline solution for about 5 minutes. Ensure that the BRAXON® fast matrix is completely covered by saline solution.
6. BRAXON® fast must be stored in a sterile solution for no more than 4 hours before implantation.
7. Ensure that the size of the BRAXON® fast matrix is suitable for containing the entire implant (see size chart).
8. Spread BRAXON® fast out on a sterile surface and place the base of the implant on the BRAXON®

fast flat disc.

9. Fold the BRAXON® fast shell on the implant so as to allow the dome-shaped part of BRAXON® fast to adapt to the roundness of the silicone implant. At this point, we recommend taking the upper part of BRAXON® fast, where the two flaps are sutured together, and lift it up so as to remove any air bubbles. This will also allow the implant to settle correctly in position.
10. Use resorbable sutures to join the anterior and posterior flaps of BRAXON® fast, following the outline of the implant. Cut out any excess part of the shell keeping enough margin ( $\geq 5$  mm) to secure the device to the pectoralis major.
11. The BRAXON® fast matrix is now ready to be implanted over the pectoralis major muscle. Use multiple sutures at the cardinal points to ensure primary stability and close contact between the device and vital, well-vascularised tissue. We recommend capitonnage sutures in the anterior portion of the device to ensure contact between the BRAXON® fast matrix and the subcutaneous tissue and prevent any gaps.

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